



BMA BIOMEDICALS



Peninsula Laboratories

Rabbit anti Nicastrin (1-17) amide

Nicastrin (abbreviated NCT) is a protein that is part of the gamma secretase protein complex, which is one of the proteases involved in processing amyloid precursor protein (APP) to the short Alzheimer's disease-associated peptide amyloid beta. Nicastrin itself is not catalytically active, but instead promotes the maturation and proper trafficking of the other proteins in the complex, all of which undergo significant post-translational modification before becoming active in the cell. It has also been identified as a regulator of neprilysin, an enzyme involved in the degradation of amyloid beta fragment.

This antibody was generated by immunization of rabbits with Nicastrin (1-17) amide coupled to a carrier protein.

TECHNICAL AND ANALYTICAL CHARACTERISTICS

Catalog Number:	T-4787 (Lot n/a)
Clonality:	Polyclonal
Host species:	Rabbit
Reactivity:	Human
Quantity:	50 µl
Format:	Neat undiluted antiserum, lyophilized, packaged under nitrogen. Reconstitute by adding 50 µl distilled water. This will give the equivalent of undiluted antiserum.
Stability and Storage:	Original vial: 1 year at 4° - 8°C from date of purchase. Stock solution and aliquots thereof: 1 year at -20°C. Avoid freeze-thaw cycles.
Immunogen:	Synthetic peptide H-Trp-Tyr-Lys-Pro-Ala-Ala-Gly-His-Ser-Ser-Tyr-Ser-Val-Gly-Arg-Ala-Ala-Gly-Leu-Leu-Ser-Gly-Leu-OH coupled to carrier protein.
Protein database:	n/a

Applications:

Tested in immunohistochemistry (IHC) and ELISA against Nicastrin (1-17) amide. Other applications like FACS or Western Blot may work as well.

Optimal dilutions should be determined by the end user.

Please see **bma.ch** for protocols and general information.

Approximate working dilution for IHC:

Frozen sections: not tested yet.

Paraffin sections: 1:500; recommended Proteinase K pre-treatment.

Related Products:

T-4790: Rabbit anti Nicastrin (656-676) (hu), neat antiserum

This product contains no preservative and is intended for laboratory use and research purposes only. Not for diagnostic or therapeutic applications.