



BMA BIOMEDICALS



Peninsula Laboratories

Y-1055 ELISA buffer Concentrate

Reconstitution Information:

Each bottle of the ELISA Buffer Concentrate contains 50ml of 20x concentrated ELISA buffer. Dilute the ELISA Buffer Concentrate to 1,000ml (1 liter) with deionized or USP pure water and mix well.

Composition

The Buffer Concentrate contains Trizma base and Trizma hydrochloride (pH 7.4), NaCl, Tween 20, 0.05% ProClin® 150 and 1% BSA.

Recommended Use

The ELISA Buffer is used for several steps in our ELISA (Enzyme Immunoassay) protocols.

Storage

Store refrigerated at 4 - 8°C.