



BMA BIOMEDICALS



Peninsula Laboratories

T-4577

Normal Rabbit Serum, pretitrated

Reconstitution Information:

Rehydrate the lyophilized neat serum with 50ml of buffer, this should be the same buffer used in the assays. It is sufficient for 500 tests.

Stability:

At least one year at 4° - 8°C from date of delivery. Minimize repeated thawing and freezing of the serum by freezing aliquots at - 20°C or below.

Related Products: T-4576: Goat anti Rabbit IgG Serum (pretitrated)
Y-1000: SEP-Columns (200mg, 50 per box)
Y-1010: SEP-Column Adapters (15 per box)
Y-1040: Buffer A for peptide extraction
Y-1045: Buffer B for peptide extraction
Y-1050: RIA Buffer

This product contains no preservative and is intended for laboratory use and research purposes only. Purchase of this product does not include authorization to use it in diagnostic or therapeutic applications.