



BMA BIOMEDICALS



Peninsula Laboratories

S-1155

Calcitonin (salmon I) ELISA [UniProt: P01263]

Calcitonin is a 32 amino acid peptide hormone that defines the calcitonin-like protein family. It is formed by the proteolytic cleavage of a larger prepropeptide, which is the product of the CALC1 gene and secreted by parafollicular cells of the thyroid gland in humans, and in the ultimopharyngeal body in animals. The hormone participates in calcium and phosphorus metabolism by reducing blood calcium, and counteracts the effects of parathyroid hormone and Vitamin D3.

This ELISA was developed with serum from rabbits immunized with Calcitonin coupled to a carrier protein.

TECHNICAL AND ANALYTICAL CHARACTERISTICS

Lot number:	A17899
Host species:	Rabbit IgG
Quantity:	96 tests
Format:	Extraction-free kit for rat serum or plasma (EIAS type).
Shelf-life:	One year from production date. Store refrigerated at 4° - 8°C.
Applications:	This ELISA has been validated with the included reagents. It is intended to be used with samples of rat origin (original protocol VII, Ab2hr.wash.Std2hr.BtON). For research use only. Please see www.bma.ch for protocols and general information.
Range:	0 - 25ng/ml
Average IC50:	0.6ng/ml
Immunogen:	Synthetic peptide H- Cys-Ser-Asn-Leu-Ser-Thr-Cys-Val-Leu-Gly-Lys-Leu-Ser-Gln-Glu-Leu-His-Lys-Leu-Gln-Thr-Tyr-Pro-Arg-Thr-Asn-Thr-Gly-Ser-Gly-Thr-Pro-NH ₂ , (Disulfide bond) coupled to carrier protein.


Cross-Reactivity:

PEPTIDE:	%:
Calcitonin (salmon I)	100
Calcitonin (eel)	100
[Asu ^{1,7}]-Calcitonin eel)	100
Calcitonin (chicken)	100
Calcitonin (human)	0
Calcitonin (rat)	0
Calcitonin C-Terminal Flanking Peptide (human)	0
Calcitonin C-Terminal Adjacent Peptide (rat)	0
Amylin (free acid) (human)	0

Cross-Reactivity (continued):

PEPTIDE:	%:
Amylin (free acid) (rat)	0
α -CGRP (human)	0
pTH Related Protein (1-34) (human)	0

Typical titration curve of Calcitonin in a competitive ELISA with this antibody:



Suggested Preparation of Standards		
	ng/ml	Range: 0.02 to 25ng/ml
Stock	1000	
S1	25.00	Add 25 μ l Stock + 975 μ l diluent
S2	6.25	Add 200 μ l S1 + 600 μ l diluent
S3	1.56	Add 200 μ l S2 + 600 μ l diluent
S4	0.39	Add 200 μ l S3 + 600 μ l diluent
S5	0.10	Add 200 μ l S4 + 600 μ l diluent
S6	0.02	Add 200 μ l S5 + 600 μ l diluent
S0	0.00	500 μ l diluent

Related Products: T-4027: anti Calcitonin (salmon I) neat antiserum
 T-4026: anti Calcitonin (salmon I) purified IgG
 T-4025: anti Calcitonin (salmon I) diluted antiserum
 S-1166: ELISA, high sensitivity for extracted samples (salmon I)

This product contains Thimerosal as a preservative and is intended for laboratory use and research purposes only. Purchase of this product does not include authorization to use it in diagnostic or therapeutic applications.

S-1155

EIAS – salmon I

30.7.2020